






OSTEOTOMES

DESIGN

THE NEW OSTEOTOMES

OTTO LEIBINGER GmbH presents a newly designed set of high quality Osteotomes developed for the lateral bone condensation of the closed sinus elevation technique. The set is consisting of straight and bayonet/angulated Osteotomes available in six diameter sizes, all with either pointed, convex or concave tips. They are also available with stop and adjusting screw. This version comes in five diameters and with convex and concave tips only. The bayonet angulated version is meant for the use in the lateral upper jaw region or in maxillary areas that are difficult to reach. Our ergonomically shaped handle fits very well into the users hand and allows good and non-exhausting working over a long period. The serration on the front part of the handle prevents from slipping while using it with gloves or wet hands. This new line of Osteotomes completes perfectly the program of instruments for implantology already successfully implemented and highly appreciated by customers all over the world.

With the same type of handle there is also a set of Chisels for bone splitting available. The Chisels come in two diameters of 4 and 7,5mm and in a straight and curved version.

	STRAIGHT [160mm	BAYONET [175mm		STRAIGHT WITH STOP [160mm	BAYONET WITH STOP [175mm
Ø 2.0mm	840/20	843/20			
Ø 2.7mm	840/27	843/27			
Ø 3.2mm	840/32	843/32			
Ø 4.2mm	840/42	843/42			
Ø 5.0mm	840/50	843/50			
Ø 6.0mm	840/60	843/60			
Ø 2.0mm	841/20	844/20		848/20	851/20
Ø 2.7mm	841/27	844/27		848/27	851/27
Ø 3.2mm	841/32	844/32		848/32	851/32
Ø 4.2mm	841/42	844/42		848/42	851/42
Ø 5.0mm	841/50	844/50		848/50	851/50
Ø 6.0mm	841/60	844/60			
Ø 2.0mm	842/20	845/20		849/20	852/20
Ø 2.7mm	842/27	845/27		849/27	852/27
Ø 3.2mm	842/32	845/32		849/32	852/32
Ø 4.2mm	842/42	845/42		849/42	852/42
Ø 5.0mm	842/50	845/50		849/50	852/50
Ø 6.0mm	842/60	845/50			



FON: +49 (0) 7463-7232

FAX: +49 (0) 7463-222

info@otto-leibinger.de

www.otto-leibinger.de